

## B1 小動物室落菌數監測結果

編號	區室	1月	2月	3月	4月	5月	6月	7月	8月	9月	10月	11月	12月
1	H20	0	0	0	0	3	0	0	0	0	1	0	0
2	H21	0	0	2	0	0	2	0	1	0	1	1	2
3	H22	0	1	0	1	0	0	0	0	0	0	1	1
4	H23	1	0	0	0	0	2	2	2	1	0	0	1
5	H24	0	0	1	0	1	2	3	1	0	2	1	1
6	H25	0	1	2	1	1	2	0	2	1	0	0	0
7	20-22 室走道	0	0	1	0	1	3	4	2	2	1	3	0
8	23-25 室走道	1	1	3	3	0	3	2	6	3	2	1	2
9	實驗室	1	1	3	1	2	0	3	5	1	1	2	0
10	洗衣室	2	1	2	2	0	5	1	2	1	3	1	1
11	滅菌室	1	0	1	0	0	0	0	2	0	0	1	0
12	清洗室	2	0	1	1	0	5	0	2	2	1	2	2
13	前走道	2	0	2	2	1	5	0	3	2	5	1	1
14	辦公室	2	0	1	1	3	2	4	4	3	4	3	3
15	更衣區	9	15	2	15	12	5	2	15	8	8	18	15
16	入口	7	4	0	15	1	8	2	2	3	16	3	3
17	20 操作台	0	0	0	0	0	0	0	0	0	0	0	0
18	21 操作台	0	0	0	0	0	0	0	0	0	0	0	0
19	22 操作台	0	0	0	0	0	0	0	0	0	0	0	0
20	23 操作台	0	0	0	0	0	0	1	0	0	0	0	0
21	24 操作台	0	0	0	1	0	0	0	0	0	0	0	0
22	25 操作台	0	1	0	0	0	0	0	0	0	0	0	1

1. BBL™ Stacker™ Plate TSB II， 5% Sheep blood (Becton, Dickinson)。
2. 方法：打開 Plate 上蓋，置於房間中央 30 分鐘，30 分鐘後蓋好上蓋，置入 37°C 培養箱中培養 48 小時，取出後計算落菌數量