

12F 動物中心各區室落菌數監測結果

編號	區室	2月	3月	5月	6月	11月	12月
1	小鼠種源室	3	3	0	0	區室整修	
2	SD繁殖室	12	12	124	160		
3	Wistar繁殖室	40	40	1	1		
4	小鼠繁殖室	3	2	0	1		
5	大鼠種源室	-	-	-	-		
6	裸鼠室	0	3	0	0		
7	儲籠室	7	5	12	10		
8	淨走廊	15	6	72	66		
9	後走廊	8	45	50	24		
10	管理室	6	4	0	0		
11	B.S.浴室	15	3	2	3		
12	潔淨1(搬運大小鼠通道)	0	28	2	9	4	4
13	潔淨2(外來加蓋繁殖)	0	0	0	2	1	1
14	潔淨3(大鼠飼育室)	2	35	29	104	216	192
15	潔淨4(大鼠飼育室)	28	58	65	76	38	36
16	潔淨5(大鼠飼育室)	0	23	9	1	6	2
17	潔淨6(小鼠飼育室)	1	3	7	7	18	11
18	潔淨區前室	0	4	4	4	29	17
19	潔淨區走廊	0	5	1	8	50	16
20	外科手術區	6	0	3	10	8	2
21	小鼠代養室1(H1國小加蓋繁殖1)	0	11	2	10	9	0
22	小鼠代養室2(H2國小加蓋繁殖2)	0	3	1	2	8	0
23	小鼠代養室H1,2前室	2	3	0	0	41	2
24	小鼠代養室3(H3)	44	97	180	91	146	34
25	大鼠代養室4(H4)	11	8	56	21	72	70
26	大、小鼠代養室H3,4前室	17	0	4	9	29	17
27	H1~H4前走道	37	10	13	14	17	4
28	基因轉殖動物代養室(H5)	2	0	2	3	16	1
29	基因轉殖動物代養室前室(H5前室)	4	1	0	5	3	10
30	感染室(H6)	0	4	4	4	1	0
31	感染室前室(H6前室)	3	12	8	5	2	3
32	H5~H6前走道	8	5	2	15	24	5
33	基因轉殖加蓋室(H7)	1	0	2	0	1	1
34	代養室(H8)	10	57	63	37	21	22
35	外來動物隔離區1(H8-1國大鼠)	1	2	3	20	14	0
36	外來動物隔離區4(H8-2國小鼠)	0	3	3	3	2	1
37	外來動物隔離區3(H8-3外來)	7	1	4	3	1	0
38	外來動物隔離區2(H8-4操作台)	0	4	2	1	0	1
39	外來動物隔離區前室(H8前室)	3	4	10	15	8	5
40	大鼠代養室H9	19	15	77	23	28	28
41	大鼠代養室H9前室	1	77	5	9	8	3
42	小鼠代養室H10	1	70	112	176	48	42
43	小鼠代養室10前室	37	3	9	5	2	13
44	代養室H8,9,10前走道	8	7	7	24	12	6
45	豬室(H11)	999	999	999	999	888	999
46	兔室1(H12)	70	24	120	94	82	49
47	豬室(H13)	999	999	999	999	777	789
48	兔室2(H14)	26	14	87	196	40	27
49	清洗室	5	7	1	10	1	9

B1 小動物室各區室落菌數監測結果

區室	1月	2月	3月	4月	5月	6月	7月	8月	9月	10月	11月	12月
H20	0	0	0	0	0	0	0	0	0	1	0	0
H21	0	1	0	0	0	0	0	0	0	0	0	0
H22	1	0	1	0	3	0	0	0	0	0	0	0
H23	0	4	0	1	1	2	1	0	0	0	0	1
H24	1	1	0	2	0	0	0	1	0	1	0	0
H25	0	1	0	0	1	0	1	0	0	0	0	0
H20-22 走道	0	0	0	0	1	0	0	0	0	0	0	0
H23-25 走道	1	0	0	0	0	0	0	0	0	0	0	0
實驗室	0	0	0	0	0	0	0	0	1	0	0	0
洗衣室	1	1	0	1	1	1	2	0	5	1	3	1
滅菌室	0	0	0	0	0	0	0	0	0	0	0	3
清洗室	0	2	1	2	1	0	0	0	0	2	0	1
前走道	1	0	4	3	2	2	1	0	3	0	1	1
辦公室	0	0	0	1	3	1	2	0	6	3	6	4
更衣室	2	3	1	2	4	11	4	5	1	1	5	10
入口	2	3	1	1	2	8	5	2	2	5	1	5
H20 操作台	0	0	0	0	0	0	0	0	0	0	0	0
H21 操作台	0	0	0	0	0	0	0	0	0	0	0	0
H22 操作台	0	0	0	0	1	0	0	0	0	0	0	0
H23 操作台	0	0	0	0	0	1	0	1	1	0	0	0
H24 操作台	1	1	0	0	0	0	1	0	0	0	0	0
25 操作台	0	0	0	0	0	0	0	0	0	0	0	1

方法：打開 Plate 上蓋，置於房間中央 30 分鐘，30 分鐘後蓋好上蓋，置入 37°C 培養箱中培養 48 小時，取出後計算落菌數量。

